





Magnetic nano-particles with no surface treatment in aqueous medinm

MAGNETIC NANO-PARTICLES DEVELOPER KIT for Biomedical application



Ferrotec Corporation (JAPAN)

1-4 Midoridaira Yokaichiba Chiba, 289-2131 JAPAN TEL 81-479-73-6752 FAX 81-479-73-6602 ffsales@ferrotec.co.jp

Ferrotec (USA) Corporation

40 Simon Street Nashua, NH 03060 USA TEL 1-603-883-9800 FAX 1-603-883-2308 sales@ferrotec.com


<u>CAUTION</u>

For Investigational Use Only. The performance characteristics of this product have not been established. This product is not sold for use in products for which prolonged contact with skin or implantation in the human body. Ferrotec Corporation and its group companies do not recommend this material as safe and effective for such uses and assumes no liability for any such use. Ferrotec Corporation and its group companies also do not warrant and specifically disclaims that the use of its ferrofluid & magnetic particles for biomedical applications covered by another party's patent do not or will not infringe the intellectual property rights of the owner of such biomedical patent.


Issued on Nov. 28, 2005 TN-EMGIIII rev.A

EMG1111 is water based slurry containing magnetic nano particles of iron oxide coated with no surfactant. The particles have a nominal diameter of about 10nm having single domain & superparamagnetic property. Therefore no hysterysis on magnetization curve can be seen as a typical data obtained by VSM (Vibrating Sampling Magnetometer). The particles also have magnetic permeability as in table and an initial susceptibility of about 0.6 typically.

The particles have about 10 nm diameter on average, however some distribution of the size can be seen as a picture of TEM (Transmittance Electron Microscope). A core of the particles are iron oxides and these are well known as compatible with living body.



VSM data for typical EMG1111 (iron oxide particles in slurry)



XRD analysis data for typical EMG1111



TEM picture for typical EMG1111

Physical properties for EMG1111 (specification or typical data)

Арреагапсе	Water based slurry
Saturation magnetization of iron oxide in water based slurry (Typical data)	About 65 emu/g
Average particle size (Typical data)	About 10 nm
Density of water based slurry at 25C	1.13~1.23 g/cm ³

From the X-ray analysis data by using XRD (X-Ray Diffraction spectroscopy), we can see the iron oxides are the mixture of Fe3O4 and gamma-Fe2O3.

EMG1111 allows the researcher to coat the particles with his/her own surfactant. The uncoated particles are also used as reference sample for some researchers.

We recommend an ultra sound treatment (sonication) prior to use of the product.

Please feel free to contact Ferrotec if you need technical assistance for the particles.